

Oxy-Care / Ute-Fix

Decubitus Prophylaxis

Oxy-Care / Ute-Fix are atraumatic and comfortable paddings for all kind of tubes and drainages that may entail skin lesions.



Oxy-Care is available in one size:

Oxy-Care REF 40 - 10 Ø 1,0 cm 7,0 cm padding length

In particular this padding protects the skin against pressure marks and irritation caused by nasal cannulae. In the range of respiration support with neonates, babies and children Oxy-Care is used as a tube padding to protect the sensitive skin of the cheeks.

Ute-Fix is available in three sizes:

Ute-Fix 1 (mini) REF 40 - 06 Ø 1,5 cm 2,0 cm padding length

Ute-Fix 2 (small) REF 40 - 05 Ø 2,0 cm 2,0 cm padding length

Ute-Fix 3 (large) REF 40 - 01 Ø 2,5 cm 2,0 cm padding length

Novo-San REF 20 - 80 120 wipes à dispenser box

This padding is especially used with tubes and drainages to reduce the risk of eudermic infection, ulceration, maceration and decubitus.

Desinfectant cleaning of tubes with Novo-San tintured wipes before attaching the adhesive surface in case of first application as well as changing the fixation assures optimal contact power.

ADVANTAGES

	patient	nursing staff	finance controller
safe, quick and comfortable handling	●	●	●
bacterial neutral – foam reduces bacteria growth	●	●	●
latex-free – hypoallergenic, protection against decubitus ulcer	●	●	●
reduction of local pressure due to low density	●		
can easily be trimmed to fit – sizes suitable for all patients	●	●	●
user-friendly – also applicable wearing gloves		●	
adapted for the use in clinics, care and Home Care	●	●	●
sterilisation is possible	●	●	●
Made in Germany	●	●	●

Oxy-Care / Ute-Fix




Oxy-Care used in Home Care

This padding protects the skin around the ear against pressure marks and irritation caused by nasal cannulae.



Oxy-Care



Pressure mark due to the tube

Ute-Fix 3 (large)

In prophylactic use it protects the skin against pressure marks or excoriation using suprapubic and transurethral catheters as well as sundry tube systems.



Ute-Fix 2 (small)

Open the padding and take the cover of the adhesive surface inside off.



Press the padding around the tube very close to the entrance of the nose. The tube has to be well cleaned. For cleaning please use e.g. **Novo-San**.



Place the tube with the padding at the entrance of the nose.