

TuBo-Grip non-adhesive

Two-Part Endotracheal Tube Holder

TuBo-Grip is a two-part endotracheal tube holder for adults, adolescent patients and pupils. The silicone layer of the elastic strap clings to the tube and offers stability even with intensely salivating patients. The fixation is softly padded in the neck and therefore specially convenient for sensitive skin.



TuBo-Grip consists of a neck padding and a middle part. The neck padding is available in two different sizes:

TuBo-Grip complete	REF 20 - 50	neck padding + middle part
TuBo-Grip complete extra long	REF 20 - 51	neck padding extra long + middle part
TuBo-Grip neck padding	REF 20 - 52	length 32 cm, width 2,5 cm
TuBo-Grip neck padding extra long	REF 20 - 53	length 47 cm, width 2,5 cm
TuBo-Grip middle part	REF 20 - 55	length 33 cm, width 1,8 cm
Novo-San	REF 20 - 80	120 wipes à dispenser

Disinfectant cleaning of the tube with **Novo-San** in case of changing the fixation assures optimal contact power of the silicone layer.

Additionally, **TuBo-Grip** is available in a set with the bite guard **TuBo-Guard**. Just contact us.

ADVANTAGES

	patient	nursing staff	finance controller	
safe, quick and comfortable handling	●	●		
safe – outstanding high stability by silicone layer + tension of elastic strap	●	●		
non-adhesive – ideal for patients with beards and sensitive skin	●	●	●	
latex-free – hypoallergenic, bacterial neutral, protection against decubitus ulcer	●	●	●	
dimensionally stable – ideal for long-term use by wide, comfortable padding	●	●		
slightly absorptive – no loss of function, readjustable at any time	●	●		☒ = do not reuse
can be cut individually – fixation fits all patients	●	●	●	☒ = latex-free
suitable also for double lumen tubes and larynx-masks		●	●	☀ = keep away from sunlight
sterilisation is possible e.g. by vapour at 134 °C	●	●	●	☂ = keep dry
economic – replace middle part, padding might stay for a longer period	●	●	●	☒ = consult instruction for use
Made in Germany	●	●	●	CE Made in Germany

TuBo-Grip

non-adhesive

NOVO

Novo Klinik-Service GmI



Metra helps helping!

Instruction for Use



When changing the fixation: Please clean the tube with **Novo-San**.

Place the neck padding behind the patient's neck, so that one end is placed over the ear and one beneath the other ear.



Place one end of the elastic strap with the silicone layer from above on the tube.



Wind the elastic strap under tension once around the tube.



Wind the elastic strap under tension once again around the tube, so that the elastic strap is crossing on the bottom side of the tube.



The optimal tension is achieved when the elastic strap is completely surrounding the tube. The soft padding lies with the white side on the skin respect. in the corner of the mouth.



Fasten the ends of the middle part under tension on the neck padding using the Velcro® strips.

Fixation in final position.